ROMÂNĂ

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Before collection

Stool sample collection kit

This kit contains a tube with a cotton swab. For an accurate test result, it is essential to follow the collection instructions below.

01 Do not open the tube until you are ready to collect the sample.

02 Do not take antibiotics, antifungal drugs, or immunosuppressive medication for **6 weeks** before sample collection.

03 If you have received radiotherapy or chemotherapy, wait at least **6 weeks** before collecting the sample.

04 Do not take probiotics for **3 weeks** before sample collection.

05 If you have undergone fertility treatment, wait **4 weeks** before collecting the sample.

06 Do not collect the sample if you have active bleeding due to hemorrhoids or menstruation, or during a flare-up of irritable bowel syndrome.

07 Maintain your usual daily habits for **2 weeks** before collection so the test reflects your “normal” microbiome.

08 Collect a sample from a normal bowel movement.

09 Keep the sample at **room temperature** after collection.

Important: do not stop prescribed medication without medical advice. Confirm collection timing with the laboratory.

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How to collect the sample

Follow each step in order.

01 Attach the barcode label to the collection tube **vertically**.

02 After a spontaneous bowel movement, wipe with a piece of toilet paper.

03 Open the tube by twisting and pulling. Make sure the swab **does not touch any surface** other than the sample.

04 Gently rub the cotton tip of the swab over the stool trace on the toilet paper. **Do not overload** the swab with stool or solid pieces; a small trace is enough for analysis.*

* Too much or too little sample material may affect or even prevent the test.

AMOUNT OF MATERIAL ON THE SWAB

✕ Too little ✓ Correct ✕ Too much

05 Immediately place the cotton end of the swab back into the tube without touching any other surface, then press to close until it **clicks**.

06 Keep the sample at **room temperature** until it is handed over to the laboratory.